



GE MRI

MRI COIL LOADING (PATIENT EFFECT) - UNDERSTANDING THE IMPACT


The patient becomes part of the RF system and changes the coil's impedance and tuning.

Good Loading → Tune & Match Easily → High SNR & Uniform Images | Poor Loading → Detune & Mismatch → Signal Loss & Artifacts

1. WHAT IS COIL LOADING?

When a patient is placed in the coil, their body (high water content) interacts with the RF magnetic field. This causes:

- Resonant frequency shift (f_0)
- Change in impedance (Z)
- Change in Q (quality factor)
- Change in SNR & penetration



The patient is part of the circuit!
Always tune & match with the patient in the coil whenever possible.


2. HOW LOADING CHANGES THE COIL

Parameter	Without Patient (Unloaded)	With Patient (Loaded)	Effect
Resonant Frequency (f_0)	Higher	Lower	↓
Impedance (Z)	Higher	Lower	↓
Quality Factor (Q)	Higher	Lower	↓
Bandwidth (BW = f_0 / Q)	Narrower	Wider	↑
SNR	Lower	Higher	↑
Penetration	Lower	Higher	↑

Typical frequency shift with body coil:
-1 to -3 MHz (on 64-128 MHz systems)


3. EQUIVALENT CIRCUIT VIEW

UNLOADED COIL
(No Patient)



High f_0 , High Z, High Q
Narrow BW

LOADED COIL
(With Patient)



Lower f_0 , Lower Z, Lower Q
Wider BW (more damping)

R_{body} represents the patient (lossy tissue) that absorbs energy and lowers Q and impedance.

4. LOADING EFFECT ON SNR & PENETRATION

GOOD LOADING
(Well matched & tuned)

- ✓ Efficient power transfer
- ✓ High SNR
- ✓ Better penetration
- ✓ Uniform images



POOR LOADING
(Detuned / Mismatched)

- ✗ Poor power transfer
- ✗ Low SNR
- ✗ Shallow penetration
- ✗ Non-uniform / noisy images



5. LOADING EXAMPLES - GE COILS

GE HD Body Coil
(Integrated)



GE 8ch Head Coil
(Receive Only)



GE Spine Array
(16 Channel)



GE Knee Coil
(8 Channel)



Different coils load differently depending on size, elements, and patient size/position.

6. TYPICAL LOADING EFFECTS (VALUES VARY BY SYSTEM & COIL)

Coil Type	Frequency Shift (Δf)	Impedance Change (ΔZ)	Q Change	Comments
Body Coil (Integrated)	-1 to -3 MHz	↓ 30-60%	↓ 2-5x	Large loading, more damping
Head Coil (8-32 ch)	-0.5 to -2 MHz	↓ 20-50%	↓ 2-4x	Moderate loading
Spine Coil (Array)	-0.5 to -1.5 MHz	↓ 20-40%	↓ 1.5-3x	Depends on FOV & patient size
Extremity Coil (Knee/Shoulder)	-0.2 to -1 MHz	↓ 10-30%	↓ 1.2-2x	Less loading

↑ = Increase ↓ = Decrease

7. LOADING VS PATIENT SIZE / POSITION

LARGER PATIENT / MORE IN COIL

- More tissue in the RF field
- More loading (lower f_0 , lower Z, lower Q)
- Higher SNR but more absorption



SMALLER PATIENT / OFF-CENTER

- Less tissue in the RF field
- Less loading (f_0 closer to unloaded)
- Lower SNR, possible hot spots



Always center the patient and use proper coil selection.
Correct positioning = consistent loading = consistent image quality.


8. HOW TO MANAGE LOADING

- ✓ Tune with patient in the coil.
- ✓ Match with patient in the coil.
- ✓ Use correct coil for exam & patient size.
- ✓ Center patient in coil (left-right, up-down).
- ✓ Use proper FOV & prescription.
- ✓ Re-tune if patient position changes.
- ✓ Use pre-scan normalization (Auto PNS, etc.).

9. SIGNS OF POOR LOADING

- ✗ Large frequency shift (> 2-3 MHz)
- ✗ Low impedance (< 10 Ω at 50 Ω system)
- ✗ Very low Q (< 10)
- ✗ Reflected power high (> -10 dB)
- ✗ Low SNR / non-uniform images
- ✗ Bands / shading / folds in image

Example: Poor Loading



Check tune, match, position, coil, and FOV.

10. LOADING CHECK - WHAT ENGINEERS MONITOR

Frequency (MHz)	Impedance (Ω)	Q Factor (Δ)	Reflected Power (dB)
62.70	52.3	20.1	-22.4
Δf (MHz)	R (Ω)	X (Ω)	
-1.35	51.8	5.2	
Tuned	Matched	Normal	Good

Ideal Reflected Power: -20 dB or lower
(Lower is better)

11. COMMON LOADING PROBLEMS & CAUSES

Problem	Possible Cause	What to Do
Large frequency shift	Wrong coil, off-center, large patient	Re-center, use correct coil, re-tune
Low impedance	Heavy loading, detuned, coil cable issue	Re-tune, check cables, inspect coil
Low Q	Over-loaded, lossy cable, coil not tuned	Re-tune, check coil & cables
High reflected power	Poor match, detuned	Re-match, check tune
Shading / bands	Off-center, poor loading, coil element failure	Re-position, coil check, run coil test

12. BEST PRACTICES SUMMARY

- ✓ Tune with patient.
- ✓ Match with patient.
- ✓ Center patient.
- ✓ Use correct coil.
- ✓ Keep coil & connectors clean & dry.
- ✓ Check cables and element indicators.
- ✓ Document loading values if out of range.

13. QUICK CHECK FLOW

- 1 Tune ($\Delta f \sim 0$ MHz with patient)
- 2 Match (Z $\sim 50 + j0 \Omega$)
- 3 Check Reflected Power (≤ -20 dB)
- 4 Check VSWR ($\leq 1.5 : 1$)
- 5 Check Q (10-30 typical)
- 6 Verify image quality



ENGINEER NOTES

- Loading is patient, coil, position, and frequency dependent.
- Always tune & match after any patient change.
- Good loading = good images, short scans, happy patients.

ABBREVIATIONS

f_0 = Resonant Frequency	Q = Quality Factor
Z = Impedance (Ω)	BW = Bandwidth
R = Resistance (Ω)	SNR = Signal-to-Noise Ratio
X = Reactance (Ω)	FOV = Field of View



SAFETY REMINDER

- High RF power can cause tissue heating.
- Ensure patient comfort and monitoring.
- Follow all site safety procedures.



THE BOTTOM LINE

Understand loading. Manage loading.
Deliver power efficiently.
Get the best SNR, uniformity, and diagnostic confidence.