


MRI GRADIENT SIGNAL CHAIN

From Pulse Sequence to Spatial Encoding

How gradient systems work – where problems occur – and how they affect the image.

PULSE SEQUENCE → WAVEFORM GENERATION → GRADIENT CONTROLLER → GPA (X,Y,Z) → GRADIENT COIL → SPATIAL ENCODING → IMAGE

1. THE GRADIENT SYSTEM SIGNAL CHAIN – OVERVIEW



WHY GRADIENTS MATTER


- Gradients encode spatial location of the MR signal.
- They turn frequency and phase into position.
- The three axes (X, Y, Z) work together to fill k-space.
- Any error in the gradient system shows up as an image artifact.

KEY GRADIENT SPECIFICATIONS


- Amplitude (Strength): mT/m
- Slew Rate: T/m/s
- Rise Time: ms (10–90%)
- Duty Cycle / Cooling: %
- Linearity: ppm (parts per million)
- Eddy Current Performance
- Acoustic Noise / Vibration
- PNS (Peripheral Nerve Stimulation) Limits

2. HOW X, Y & Z GRADIENTS CREATE AN IMAGE

ROLES OF THE GRADIENTS



3D VIEW – ORIENTATION OF AXES




3. GRADIENT WAVEFORMS – EXAMPLES (X AXIS)

Function	Waveform (Ideal)	Purpose
Slice Select (Gz)		Select slice thickness.
Phase Encode (Gy)		Imparts phase variation.
Readout (Gx)		Encodes frequency (left → right).
Spoiler / Crusher		Dephases residual transverse magnetization.

Real waveforms are not ideal – they have finite rise time and system delays.

4. GRADIENT POWER AMPLIFIER (GPA) – BLOCK DIAGRAM



5. WHERE PROBLEMS OCCUR IN THE CHAIN

Stage	What Can Go Wrong	What You May See (Image / System)	How to Check
Waveform Generation	Incorrect scaling, timing, truncation	Wrong FOV, distortion, shift, incomplete k-space	Protocol review, logs, comparison scan
Gradient Controller	Communication error, scaling, timing jitter	Intermittent artifacts, dropouts, wrong gradients	Error logs, oscilloscope, self test
GPA (X/Y/Z)	Overcurrent, overheating, amplifier fault	System errors, no image, axis not working	Service Browser logs, GPA diagnostics
Gradient Coil	Open/short, insulation fault, loose connection	Axis offline, high noise, unusual heating	Coil resistance/inductance, insulation test
Cables / Connectors	Loose, corroded, intermittent	Intermittent artifacts, axis errors	Visual inspect, wiggle test, continuity test
System / Cooling	High temp, low flow, fans, filters	Thermal faults, reduced performance	Temp logs, flow, fans, water quality

6. GRADIENT HARDWARE PATH



Gradients are water-cooled. Proper cooling is critical for performance and longevity.

7. COMMON IMAGE ARTIFACTS FROM GRADIENT PROBLEMS

Artifact (What You See)	Likely Cause	Which Axis Usually Involved	Example
Geometric Distortion / Warping	Gradient nonlinearity, eddy currents	Usually Readout (X)	
Image Shift	Gradient offset (DC), incorrect timing	Any Axis	
N/2 Ghosting (EPI)	Polarity mismatch, timing error	Phase Encode (Y)	
Scaling Error (Wrong FOV)	Amplitude calibration error	Any Axis	
Blurring / Poor Resolution	Low amplitude, long rise time	Usually Readout (X)	
Directional Streaking	Eddy currents, poor pre-emphasis	Along the Readout Axis	

8. GRADIENT TROUBLESHOOTING APPROACH

- Review the image and recognize the pattern.
 - Check the system and gradient status for any errors.
 - Determine which axis is likely involved.
 - Change orientation (swap freq/phase) and see if the artifact follows the axis.
 - Run prescan / coil check / system baseline.
 - Check gradient calibration (offset, gain, linearity, delay).
 - Use test tools (impulse, ramp, PNS, dB/dt monitor).
 - Isolate the stage (controller → GPA → coil → cables).
- GOOD IMAGES REQUIRE**
- Correct waveforms
 - Accurate current
 - Fast response
 - Linear fields
 - Proper timing
 - Effective eddy current compensation
 - Stable temperature and cooling
- When gradients are accurate in time, amplitude and linearity – the image is sharp, undistorted and reliable.



ENGINEER NOTES

- Always document changes and results.
- Compare with baseline and other coils/axis.
- Escalate if data is inconclusive.

ABBREVIATIONS

GPA = Gradient Power Amplifier	RO = Readout	dB/dt = Rate of Change of Gradient Field
FOV = Field Of View	PNS = Peripheral Nerve Stimulation	EPI = Echo Planar Imaging
PE = Phase Encode	DC = Direct Current	ADC = Analog to Digital Converter



SAFETY REMINDER

- High currents and voltages are present in GPA.
- Allow time for discharge. Follow lockout/tagout.
- Verify cooling. Hot gradients can fail quickly.
- Follow all GE safety procedures.



THE BOTTOM LINE

Gradients encode space.
Good gradients make great images.
Understand the chain. Verify the system. Fix the cause.



GE HealthCare
Imagination at work