

WHEN THE TECHNOLOGIST SAYS "THE IMAGES DON'T LOOK RIGHT..."

A SYSTEMATIC APPROACH FOR MRI ENGINEERS



Use this step-by-step roadmap to investigate image quality issues from the applications perspective. Start simple → end complex.

1		PROTOCOL & CLINICAL GOAL - Is the correct exam and sequence prescribed? - Are the right views/planes included? - Is contrast required? Was it given? - Review prior studies for comparison.	Ask: What are we trying to show? Are we looking at the right thing?
2		PARAMETERS & SCAN SETTINGS - TR / TE / TI appropriate for the goal? - FOV, matrix, slice thickness, and bandwidth? - Fat sat or contrast timing correct? - Any recent protocol or parameter changes?	Ask: Are the parameters optimized for the clinical goal?
3		PATIENT POSITIONING & IMMOBILIZATION - Is the patient centered and aligned? - Are landmarks and prescription planes correct? - Could motion (voluntary/involuntary) be the cause? - Is breath-hold or gating working properly?	Ask: Is the patient in the right place and holding still?
4		COIL & HARDWARE CONNECTIONS - Is the correct coil selected and connected? - Any coil element failure (noise, stripes, dropouts)? - Coil is centered and snug (no gaps or tilt)? - Check cable connections and locking.	Ask: Is the right coil working as it should?
5		PRESCAN / CALIBRATION / CENTRING & FREQUENCY (CF) - Were prescan/calibrations completed successfully? - Is the FOV properly centered on the anatomy? - Is the frequency correct (water peak)? - Any aliasing or wrap-related problems?	Ask: Is the system measuring the right signal?
6		SHIM & B₀ HOMOGENEITY - Poor uniformity can cause distortion, signal loss. - Review auto shim results and map. - Any metal, air gaps, or patient size issues?	Ask: Is the magnetic field as uniform as possible?
7		RF SYSTEM (B₁) & POWER - Is RF gain appropriate (not too low or too high)? - Any dielectric shading, poor fat sat, or inhomogeneity? - Check SAR limits, duty cycle, and coil tuning/match.	Ask: Is the RF doing its job consistently?
8		GRADIENT SYSTEM - Any streaking, geometric distortion, or ghosting? - Check for gradient performance or heating issues. - Review eddy current compensation if applicable.	Ask: Are the gradients performing as expected?
9		SYSTEM, SOFTWARE & WORKFLOW - Software version, known service bulletins? - Hardware faults or error logs? - Network / archive / workstation issues?	Ask: Is everything else in the system healthy?

QUICK TRIAGE: MATCH THE PROBLEM TO LIKELY CAUSES

What You See	Possible Causes	Start Here
 Blurry / Low Detail	Motion, wrong FOV/matrix, low SNR, thick slices	2 3 5
 Signal Loss (Dark)	Wrong TI/TE, off-resonance, B ₀ inhomogeneity, poor coil	2 6 4
 Too Bright / Saturated	Saturation (flow/in-plane), TR too short, gain too high	2 5 7
 Ghosting	Motion, flow, wrap/aliasing, phase encode problem	3 5 8
 Streaking / Distortion	Metal, B ₀ issues, gradient nonlinearity, eddy currents	6 8 2
 Incomplete Fat Sat	B ₁ inhomogeneity, wrong method, off-center	7 4 5
 Venous Enhancement	Contrast timing, TR too short, flow vs tissue timing	5 2 1
 Poor MRA Arteries	Saturation, slow flow, timing, suboptimal method	2 6 1

ARTIFACTS AT A GLANCE

Motion  - Blurry, double edges - Ghosting - Phase/Wrap Check: 3, 5	Metal Susceptibility  - Signal void - Distortion - Streaking Check: 6, 8
Chemical Shift  - Bright/dark edges - More at high BW - Opposite sides Check: 2	Truncation (Wrap)  - Anatomy appears on opposite side - FOV too small Check: 2, 5
Flow / In-Plane Sat  - Loss of flow signal - Time-of-Flight sensitive Check: 2, 5, 1	Dielectric Shading (B₁)  - Uneven fat sat - Signal variation - Peripheral loss Check: 7
Gibbs (Truncation)  - Ripples near sharp edges - High contrast boundaries Check: 2	Nyquist / Aliasing  - Repeated anatomy - FOV too small Check: 2, 5

GADOLINIUM CHECKLIST

- ☒ Correct patient (name, weight, allergies, eGFR)?
- ☒ Correct dose and concentration?
- ☒ IV access patent; power injector working?
- ☒ Saline flush before and after injection?
- ☒ Timing/delay correct for sequence?
- ☒ Trigger (bolus/auto-trigger) placed correctly?
- ☒ Any adverse reaction or extravasation?
- ☒ Document lot #, dose, time, and site.



Tips

- For dynamic studies, consistency is key.
- Review prior images for timing.
- Arterial phase = 15–25 sec (brain/neck)
20–30 sec (body).



REMEMBER

- Start simple, think systematically.
- Most image problems are applications or setup issues, not hardware failures.
- The best engineers solve problems before the patient leaves.

SNR QUICK WINS



Increase NEX / averages



Increase voxel size (thickness or matrix)



Increase bandwidth



Use phased array coils



Optimize parameter match for the goal



Reduce motion

GOLDEN RULES

- ★ Quality in → Quality out.
- ★ Garbage in → Garbage out.
- ★ Change one thing at a time.
- ★ Verify before you escalate.
- ★ Document what you did.
- ★ Communicate with respect.
- ★ You are the engineer—be the problem solver.



KEY TAKEAWAY

Follow the roadmap. Ask the right questions.
Check the basics first. Solve problems together.
Great images = Great answers = Great patient care.