



1.15 CHAPTER 1 SUMMARY – FROM ATOM TO MR IMAGE

Putting it all together: How MRI signal is created, evolves, and becomes an image.

THE MRI SIGNAL FORMATION JOURNEY – STEP BY STEP



TIME DOMAIN SIGNALS			
	FID (After 90° Pulse)	SPIN ECHO (Refocused)	GRADIENT ECHO (GRE) (Gradient Reversed)
Rephasing Mechanism	None	180° RF pulse	Gradient reversal
Decay Constant	T_2^* (includes ΔB_0)	T_2 (no ΔB_0)	T_2^* (includes ΔB_0)
Echo Time (TE)	—	At TE	At TE
Signal Amplitude	Highest at $t = 0$	Higher (longer TE)	Lower (faster decay)
Sensitivity to Field Inhomogeneity	High	Low (static only removed)	High
Typical Use	Not used for imaging	SE, TSE, FSE	GRE, FLASH, SPGR



SEQUENCE COMPARISON AT A GLANCE		
Feature	Spin Echo (SE)	Gradient Echo (GRE)
Refocusing Method	180° RF pulse	Gradient reversal
180° Pulse	Yes	No
Echo Type	Spin Echo	Gradient Echo
Decay Constant	T_2	T_2^*
Sensitivity to ΔB_0 (Inhomogeneity)	Low	High
SNR (for same TE)	Higher	Lower
Imaging Speed	Slower	Faster
Tissue Contrast	T_2 -weighted (via TE)	T_2^* -weighted (via TE)
Common Sequences	SE, TSE, FSE	GRE, FLASH, SPGR, MPRAGE
Best For	High SNR, robust to inhomogeneity	Fast imaging, short TE, susceptibility contrast (e.g., SWI)

KEY TAKEAWAYS	
✓	MRI begins with the hydrogen proton and its spin in B_0 .
✓	RF excitation tips magnetization into the transverse plane.
✓	Precessing spins induce a voltage in the receive coil.
✓	Spins dephase due to T_2 (true) and T_2^* (field inhomogeneity).
✓	FID is the raw signal right after excitation.
✓	180° pulse refocuses static field inhomogeneity $\rightarrow$ Spin Echo.
✓	Gradient reversal refocuses spins $\rightarrow$ Gradient Echo.
✓	Echo amplitude decays: SE with T_2 , GRE with T_2^* .
✓	Spatial encoding with gradients fills k-space.
✓	Fourier Transform converts k-space data to an image.
✓	Contrast depends on sequence parameters and tissue properties (T_1 , T_2 , PD, T_2^*).

ENGINEERING INSIGHT

- Shim quality (B_0 uniformity) directly affects T_2^* and GRE.
- Gradient performance (linearity, timing) is critical for k-space accuracy and image quality.
- RF system (amplitude, phase, stability) determines flip angle accuracy and signal consistency.
- Coil performance (preamp, tuning, decoupling, Q) affects SNR and image uniformity.
- Environmental factors (temperature, vibrations) can impact stability and image quality.

FIELD SERVICE TIPS

- Verify shim before advanced sequences.
- Check gradient calibration and timing.
- Confirm RF output power and flip angle calibration.
- Inspect RF cables, preamps, and coil connections.
- Monitor for eddy current or vibration issues.
- Use phantom scans to baseline performance.
- Document and trend key QA metrics.

REAL-WORLD ANALOGY

Think of a group of runners starting together.

No 180° Pulse (Raw FID)



They spread out and never return.

Spin Echo (180° Pulse)



The 180° pulse makes them turn around and meet again at TE (echo).

Gradient Echo (GRE)



Changing the path (gradient) makes them meet at TE.

WHY IT MATTERS

- SE gives reliable images with high SNR.
- GRE enables fast imaging and specialty contrasts (susceptibility, flow, functional).
- Understanding this physics helps you troubleshoot, optimize, and maintain every MRI system.
- Good physics knowledge = better service, better uptime, better patient care.

QUICK REFERENCE

$\gamma/2\pi$ (proton) 42.576 MHz/T
Larmor Frequency $f_0 = \gamma/2\pi \cdot B_0$
 T_1 = longitudinal relaxation time
 T_2 = transverse relaxation time (spin-spin)
 T_2^* = apparent transverse relaxation time (includes ΔB_0)
TE = time to echo
k-space = spatial frequency domain
Image = 2D Fourier Transform of k-space data

